



Dual Panel Digital Upgrade

DP4/4+

Digital Upgrade Kit for AMX 4/4+

Powered by Canon DR Technology



Maximize
Mobile Productivity



Digital Radiography Beyond Traditional Applications

Expand the purpose of your AMX mobile with a Dual Panel upgrade! Continue to tackle general bedside exams using the 7.6 lbs Canon CXDI-55 Series digital flat panel with a 14"x17" imaging area, then compliment your upgrade with a Canon CXDI-60 Series flat panel that boasts a mere 5.9 lbs and 9"x 11" imaging area. Ideal for NICU and pediatrics, the CXDI-60C provides less X-ray dosage to the patient thanks to its Cesium (CsI) scintillator.

FEATURES / SPECIFICATIONS

- Create, manage, manipulate digital images at bedside
- No more carrying cassettes to a CR reader
- Flexible patient positioning
- Improved workflow and patient throughput
- An efficient, inexpensive alternative
- Minimal downtime / upgrade in hours / 2nd day training
- 3-5 second image preview
- Full DICOM connectivity with storage, print class and modality work list with MPPS support

GENERATOR

Compatible Models	AMX 4, AMX 4+
kW Rating	12.5
kVp Range	50 to 125 kVp
mA Range	100 mA Unselectable
# of kVp Stations	24
# of mA Stations	30
mAs Range	0.4 - 320 mAs
Generator Integration	No (2 PT Technique)
Battery / Capacitor	Battery Only

X-RAY TUBE

Anode Heat Capacity	140 kHu
Focal Spot	0.8 mm
Target Angle	14°
kHu / kVp	140 / 125

MOBILE UNIT

Weight	1225 lbs
Lock Release Points	2
Size (L x W x H)	45" x 25" x 76" in (114.3 x 63.5 x 193 cm)
Column Rotation	+ / - 270°
Reach	40" in (101.6 cm)
Tube Travel	Vertical: 52.5" in (133.4 cm), Horizontal: 24" (61 cm)
Motor Drive	Variable Speed Drive

OPTIONS

Digital Detector	CXDI-50 Series, CXDI-55 Series, CXDI-60 Series
Grid	<i>Please contact authorized representative</i>

Certifications: cGMP, ISO9001



VIRTUAL IMAGING, INC.

A Canon U.S.A. Company
 720 S. Powerline Rd, Suite E, Deerfield Beach, FL 33442
 PH: (800) 428-6191, (954) 428-6191 FX: (954) 428-6195
www.virtualimaging-fl.com